



COMBIFLUX Laboratory pump for low volumes

New!

One motor – several tube sets: even more universal and faster





One for all – versatile motor for different tube sets



Now the approved JUNIORFLUX pumps are also available with a separate tube set and motor. So, you only need one motor that you can mount on differing tube sets with a special quick action coupling within a few seconds. This way your work will be even more universal and you will save precious time and cost.

The COMBIFLUX pumps now offer you even more benefits for filling, emptying, transferring and dispensing small amounts of liquid out of hobbocks, canisters, narrow-necked containers and barrels.



- Special quick action coupling for the safe and fast mounting of the motor on different tube sets
- High safety factor: powerful 230 Watt motor with built-in overload cut-out switch, double insulated and splash-proof
- Easy cleaning: outer tube can be removed by hand
- 2-speed switch for sensitive filling operations
- Ergonomic handle
- Pumps sealless in Polypropylene (PP), Polyvinylidene fluoride (PVDF) and Stainless steel 1.4571 (S)
- Integral clip to suspend the pump
- Hook on pump outlet for fixing a hand tap
- With a weight of only 2 kg and being easy to handle it can be used wherever small amounts of low flammability liquids, like acids, alkaline solutions, solvents, detergents, must be decanted.



FLUX-GERÄTE GMBH Talweg 12 · D-75433 Maulbronn Tel. 07043 101-440 · Fax 07043 101-444 Fax International +49 7043 101-555 info@flux-pumpen.de · www.flux-pumpen.de