STENNER PUMPS

THE M128 SERIES PERISTALTIC METERING PUMP

Reliable injection of 1 ounce of solution to 128 ounces of water





Specifically developed to treat livestock water with:

- Antibiotics
- Vaccines
- 🗸 Dietary supplements
- 🗸 Acids
- Oxidizers (iodine, chlorine, chlorine dioxide, peracetic acid and hydrogen peroxide)
- 🗸 Disinfectants

M128 SERIES BENEFITS:


- · Accurate
- Easy to install: settings are pre-set at the factory
- Economical to maintain: simple tube replacement and one piece injection check valve
- · No need for filters upstream to pull out grit, sand and debris
- No need for back pressure valves
- Optional panel system for turnkey installation



THE M128 SERIES PERISTALTIC METERING PUMP

The M128 injects at a fixed amount per pulse which is actuated by a dry contacting water meter. The pump has a cable with four lead wires, two that must be connected to the water meter. Installations with a house controller can utilize the two remaining wires to register water consumption to the controller.

BUILD ITEM NUMBER



To build item number, select code based on application requirements. Item number is required for placing orders.

Pump Tube	05	#5 tube	
	07	#7 tube	
NOTE: Refer to sizing chart			

Voltage/Hertz	1	120V 60Hz
	2	220V 60Hz
Pulse Setting	1	1 gallon/pulse
	2	2 pulses/gallon
	3	1 liter/pulse (not available with panel system)
	4	4 pulses/gallon
	5	10 pulses/gallon (not available with 95 gpd pump)
System (optional)	Р	Pump and corresponding meter mounted on panel



M128 Pump



M128 System

SIZING FOR 1:128 DOSING

Model	Maximum Pressure	Pump Tube Number		Output ⁄60 Hz ml/min	Max. System for 1:128 gal/min	
M05	25 psi (1.7 bar)	#5	8.4	249.7	8.4	31.8
M07*	100 psi (6.9 bar)	#7	3.5	105.1	3.5	13.2

NOTE: To inject at a ratio of 1:128 at flow rates higher that 3.5 gpm up to 100 psi, the injection point can be installed after a pressure regulator. This will allow flow rates up to 8.4 gpm up to 25 psi. Contact factory for assistance with sizing

* Injection check valve included with pumps rated 26-100 psi (1.8-6.9 bar).

ALTERNATIVE SIZING FOR M07 MODEL USING 1 Pulse/Gallon METER

Dosage oz/gal	Dosage Ratio	Max. System Water Flow gal/min	Pulse Setting
1.0	1:128	3.5	1
0.50	1:256	7.0	2
0.25	1:512	14.0	4

NOTE: 1oz=29.57 ml; 1 gallon=3.78 liters; Contact factory for assistance with sizing

NOTICE: The information within these charts is solely intended for use as a guide. The output data is an approximation based on pumping water under a controlled testing environment. Many variables can affect the output of the pump. Stenner Pump Company recommends that all metering pumps undergo field calibration by means of analytical testing to confirm their outputs.

The information contained in this flyer is not intended for specific application purposes. Stenner Pump Company reserves the right to make changes to prices, products, and specifications at any time without prior notice. M128F0411



Stenner Pump Company 3174 DeSalvo Road Jacksonville, Florida 32246 USA Phone US Toll Free Fax

904.641.1666 800.683.2378 904.642.1012 www.stenner.com sales@stenner.com © Stenner Pump Company All Rights Reserved